





Product data sheet REF 09.87.55

03.06.2026

| General data | | | |
|--|---|-------------------------------------|--|
| article number | 09.87.55 | liposuct cann 15cm Ø 3.0mm hdl conn | |
| Unit | piece | Customs code (HS code) | 90189084 |
| GMDN | 45414 | GTIN | 04046826014391 |
| Risk category | 2a (EU MDD) | UMDNS-Code | 17-640 |
| Classification | 6 (EU MDD) | Base UDI-DI | |
| Net weight (kg) | 0,056 | Operating instructions | G098395 |
| Latex free | yes | Preparation instructions | MPK11 |
|  | | | |
| Medicon media | 450.00.60 (catalogue plastic surgery) 451.03.82 (brochure liposuction) | | |
| Complete picture | | Detail images | |
|  | | | |
| size | 15cm Ø 3.0mm | Coating | Artificially created galvanic gold layer (max 0.04mm) on parts of the surface of instruments |
| Storage conditions | | | |
| No special requirements | | | |
| Materials | | | |
| | Material number (DIN designation) | | |
| Main part | 1.4301 (X5CrNi18-10) | | |
| Handle | 1.4305 (X8CrNiS18-9) | | |