



Product data sheet REF 09.87.73

11.12.2025

General data			
article number	09.87.73	Entner liposuct cann ang 20cm Ø 3mm 22 holes hdl conn	
Unit	piece	Customs code (HS code)	90189084
GMDN	45414	GTIN	04046826014568
Risk category	2a (EU MDD)	UMDNS-Code	17-640
Classification	6 (EU MDD)	Base UDI-DI	
Net weight (kg)	0,057	Operating instructions	
Latex free	yes	Preparation instructions	MPK11
			
Medicon media	450.00.60 (catalogue plastic surgery) 451.03.82 (brochure liposuction)		
Complete picture		Detail images	
			
size	20cm Ø 3mm	Coating	Artificially created galvanic gold layer (max 0.04mm) on parts of the surface of instruments
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4301 (X5CrNi18-10)		
Handle	1.4305 (X8CrNiS18-9)		