



Product data sheet REF 09.87.84

25.08.2026

General data			
article number	09.87.84	transfer adaptor LL/w LL/w	
Unit	piece	Customs code (HS code)	90189084
GMDN	35384	GTIN	04046826014674
Risk category	2a (EU MDD)	UMDNS-Code	11-729
Classification	2 (EU MDD)	Base UDI-DI	
Net weight (kg)	0,004	Operating instructions	
Latex free	yes	Preparation instructions	MPK02
			
Medicon media			
Complete picture		Detail images	
size		Coating	Electroplated chrome layer on the surface
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	2.0401 (CuZn39Pb3)		