



Product data sheet REF 09.87.87

11.12.2025

General data			
article number	09.87.87	Entner liposuct cann bay 30cm Ø 4mm 30 holes hdl conn	
Unit	piece	Customs code (HS code)	90189084
GMDN	45414	GTIN	04046826014704
Risk category	2a (EU MDD)	UMDNS-Code	17-640
Classification	6 (EU MDD)	Base UDI-DI	
Net weight (kg)	0,020	Operating instructions	
Latex free	yes	Preparation instructions	MPK11
			
Medicon media	450.00.60 (catalogue plastic surgery) 451.03.82 (brochure liposuction)		
Complete picture		Detail images	
			
size	30cm Ø 4mm	Coating	Artificially created galvanic gold layer (max 0.04mm) on parts of the surface of instruments
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part Handle	1.4301 (X5CrNi18-10)		
	1.4305 (X8CrNiS18-9)		