





## Product data sheet REF 09.88.45

03.06.2026

General data			
article number	09.88.45	hdl w tube conn rotating f liposuct cann	
Unit	piece	Customs code (HS code)	90189084
GMDN	45414	GTIN	04046826014902
Risk category	2a (EU MDD)	UMDNS-Code	17-640
Classification	6 (EU MDD)	Base UDI-DI	
Net weight (kg)	0,138	Operating instructions	G098395
Latex free	yes	Preparation instructions	MPK03
			
Medicon media	450.00.60 (catalogue plastic surgery ) 451.03.82 (brochure liposuction )		
Complete picture		Detail images	
			
size		Coating	Artificially created galvanic gold layer (max 0.04mm) on parts of the surface of instruments
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4305 (X8CrNiS18-9)		