



Product data sheet REF 09.89.15

25.08.2026

General data			
article number	09.89.15	lipo injection cann Ø 0.70x40mm LL conn	
Unit	piece	Customs code (HS code)	90189084
GMDN	45414	GTIN	04046826436353
Risk category	2a (EU MDD)	UMDNS-Code	17-640
Classification	6 (EU MDD)	Base UDI-DI	
Net weight (kg)	0,004	Operating instructions	G098395
Latex free	yes	Preparation instructions	MPK16
			
Medicon media	450.00.60 (catalogue plastic surgery)		
Complete picture		Detail images	
			
size	Ø 0.70x40mm	Coating	Artificially created galvanic gold layer (max 0.04mm) on parts of the surface of instruments
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4301 (X5CrNi18-10)		
Others	1.4305 (X8CrNiS18-9)		