





Product data sheet REF 09.89.26

03.06.2026

| General data | | | |
|--|--|--|--|
| article number | 09.89.26 | von Gregory liposuct cann 18cm sh 3x180mm hole Ø 1.3mm | |
| Unit | piece | Customs code (HS code) | 90189084 |
| GMDN | 45414 | GTIN | 04046826444600 |
| Risk category | 2a (EU MDD) | UMDNS-Code | 17-640 |
| Classification | 6 (EU MDD) | Base UDI-DI | |
| Net weight (kg) | 0,018 | Operating instructions | G098395 |
| Latex free | yes | Preparation instructions | MPK11 |
|  | | | |
| Medicon media | 450.00.60 (catalogue plastic surgery) | | |
| Complete picture | | Detail images | |
|  | | | |
| size | 18cm 3x180mm Ø 1.3mm | Coating | Artificially created galvanic gold layer (max 0.04mm) on parts of the surface of instruments |
| Storage conditions | | | |
| No special requirements | | | |
| Materials | | | |
| | Material number (DIN designation) | | |
| Main part | 1.4301 (X5CrNi18-10) | | |
| Handle | 1.4301 (X5CrNi18-10) | | |
| Others | 1.4021 (X20Cr13) | | |