



Product data sheet REF 09.95.11

25.08.2026

General data			
article number	09.95.11	traction cable only stainless flex 90cm	
Unit	piece	Customs code (HS code)	90189084
GMDN	11268	GTIN	04046826014995
Risk category	1r (EU MDR)	UMDNS-Code	13-828
Classification	6 (EU MDR)	Base UDI-DI	4046826;B;R001;02;ag;37
Net weight (kg)	0,017	Operating instructions	GATDR00X
Latex free	yes	Preparation instructions	MPK01-1
			
Medicon media	450.00.03 (general catalogue the art of surgery) 450.00.53 (catalogue cardio)		
Complete picture		Detail images	
			
size	90cm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4301 (X5CrNi18-10)		
Others	1.4305 (X8CrNiS18-9)		