



Product data sheet REF 09.95.18

25.08.2026

General data			
article number	09.95.18	metal probe tip only Ø 3mm	
Unit	piece	Customs code (HS code)	90189084
GMDN	11268	GTIN	04046826015022
Risk category	1r (EU MDR)	UMDNS-Code	13-828
Classification	6 (EU MDR)	Base UDI-DI	4046826;B;R001;02;ag;37
Net weight (kg)	0,001	Operating instructions	GATDR00X
Latex free	yes	Preparation instructions	MPK01-1
			
Medicon media	450.00.03 (general catalogue the art of surgery) 450.00.53 (catalogue cardio)		
Complete picture		Detail images	
			
size	Ø 3mm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4305 (X8CrNiS18-9)		