



Product data sheet REF 09.96.54

03.06.2026

General data			
article number	09.96.54	Schmid vasc irrig cann Ø 4mm 15cm LL	
Unit	piece	Customs code (HS code)	90189084
GMDN	10730	GTIN	04046826015343
Risk category	2a (EU MDD)	UMDNS-Code	10-564
Classification	6 (EU MDD)	Base UDI-DI	
Net weight (kg)	0,016	Operating instructions	G098395
Latex free	yes	Preparation instructions	MPK02
			
Medicon media	450.00.03 (general catalogue the art of surgery) 450.00.53 (catalogue cardio)		
Complete picture		Detail images	
			
size	Ø 4mm 15cm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4301 (X5CrNi18-10)		
Others	1.4305 (X8CrNiS18-9)		