



Product data sheet REF 09.99.01

25.08.2026

General data			
article number	09.99.01	PEN suction cannula hdl	
Unit	piece	Customs code (HS code)	90189084
GMDN	38749	GTIN	04046826434229
Risk category	2a (EU MDD)	UMDNS-Code	10-566
Classification	6 (EU MDD)	Base UDI-DI	
Net weight (kg)	0,027	Operating instructions	G098395
Latex free	yes	Preparation instructions	MPK02
			
Medicon media	450.00.58 (catalogue neuro cranialis instr f cranial neurosurgery) 450.00.59 (catalogue neuro spinalis instr f spinal neurosurgery) 450.00.62 (catalogue Medicon instruments Oral and Maxillofacial Surgery) 451.02.56 (brochure PEN)		
Complete picture		Detail images	
			
size		Coating	Anodizing is the process of enriching titanium with a hard anodized layer on aluminum components.
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4301 (X5CrNi18-10)		
Handle	3.3215 (AlSi1MgMn)		
Others	1.4305 (X8CrNiS18-9)		