



## Product data sheet REF 10.00.12

03.06.2026

| General data                                                                                                                                                                  |                                                                                |                                  |                         |
|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------|----------------------------------|-------------------------|
| article number                                                                                                                                                                | 10.00.12                                                                       | Medicon Halsted nh 11.5cm smooth |                         |
| Unit                                                                                                                                                                          | piece                                                                          | Customs code (HS code)           | 90189084                |
| GMDN                                                                                                                                                                          | 12726                                                                          | GTIN                             | 04046826015473          |
| Risk category                                                                                                                                                                 | 1r (EU MDR)                                                                    | UMDNS-Code                       | 12-726                  |
| Classification                                                                                                                                                                | 6 (EU MDR)                                                                     | Base UDI-DI                      | 4046826;B;R002;01;ah;3J |
| Net weight (kg)                                                                                                                                                               | 0,020                                                                          | Operating instructions           | GATDR00X                |
| Latex free                                                                                                                                                                    | yes                                                                            | Preparation instructions         | MPK23-1                 |
| <div></div> |                                                                                |                                  |                         |
| Medicon media                                                                                                                                                                 | 450.00.03 (general catalogue the art of surgery)<br>450.00.37 (catalogue ENT ) |                                  |                         |
| Complete picture                                                                                                                                                              |                                                                                | Detail images                    |                         |
|          |                                                                                |                                  |                         |
| size                                                                                                                                                                          | 11.5cm                                                                         | Coating                          |                         |
| Storage conditions                                                                                                                                                            |                                                                                |                                  |                         |
| No special requirements                                                                                                                                                       |                                                                                |                                  |                         |
| Materials                                                                                                                                                                     |                                                                                |                                  |                         |
|                                                                                                                                                                               | Material number (DIN designation)                                              |                                  |                         |
| Main part                                                                                                                                                                     | 1.4021 (X20Cr13)                                                               |                                  |                         |