



## Product data sheet REF 10.26.73

25.08.2026

| General data                                                                                                                                                        |                                                                                                                 |                                                                                     |                                                                                              |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------|
| article number                                                                                                                                                      | 10.26.73                                                                                                        | Ryder (French Eye) HM nh 23cm                                                       |                                                                                              |
| Unit                                                                                                                                                                | piece                                                                                                           | Customs code (HS code)                                                              | 90189084                                                                                     |
| GMDN                                                                                                                                                                | 12726                                                                                                           | GTIN                                                                                | 04046826016487                                                                               |
| Risk category                                                                                                                                                       | 1r (EU MDR)                                                                                                     | UMDNS-Code                                                                          | 12-726                                                                                       |
| Classification                                                                                                                                                      | 6 (EU MDR)                                                                                                      | Base UDI-DI                                                                         | 4046826;B;R002;01;ah;3J                                                                      |
| Net weight (kg)                                                                                                                                                     | 0,077                                                                                                           | Operating instructions                                                              | GATDR00X                                                                                     |
| Latex free                                                                                                                                                          | yes                                                                                                             | Preparation instructions                                                            | MPK23-1                                                                                      |
|   |                                                                                                                 |                                                                                     |                                                                                              |
| Medicon media                                                                                                                                                       | 450.00.03 (general catalogue the art of surgery)<br>450.00.37 (catalogue ENT )<br>450.00.53 (catalogue cardio ) |                                                                                     |                                                                                              |
| Complete picture                                                                                                                                                    |                                                                                                                 | Detail images                                                                       |                                                                                              |
|                                                                                   |                                                                                                                 |  |                                                                                              |
| size                                                                                                                                                                | 23cm                                                                                                            | Coating                                                                             | Artificially created galvanic gold layer (max 0.04mm) on parts of the surface of instruments |
| Storage conditions                                                                                                                                                  |                                                                                                                 |                                                                                     |                                                                                              |
| No special requirements                                                                                                                                             |                                                                                                                 |                                                                                     |                                                                                              |
| Materials                                                                                                                                                           |                                                                                                                 |                                                                                     |                                                                                              |
|                                                                                                                                                                     | Material number (DIN designation)                                                                               |                                                                                     |                                                                                              |
| Main part                                                                                                                                                           | 1.4021 (X20Cr13)                                                                                                |                                                                                     |                                                                                              |
| Others                                                                                                                                                              | K40 (WcCo)                                                                                                      |                                                                                     |                                                                                              |