



Product data sheet REF 13.11.75

03.06.2026

General data			
article number	13.11.75	Schmieden lig needle d/e 24cm	
Unit	piece	Customs code (HS code)	90189084
GMDN	32864	GTIN	04046826019334
Risk category	1r (EU MDR)	UMDNS-Code	15-234
Classification	6 (EU MDR)	Base UDI-DI	4046826;B;R001;03;ag;3J
Net weight (kg)	0,033	Operating instructions	GATDR00X
Latex free	yes	Preparation instructions	
			
Medicon media	450.00.03 (general catalogue the art of surgery)		
Complete picture		Detail images	
  			
size	24cm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4021 (X20Cr13)		