








Product data sheet REF 13.32.05

11.12.2025

| General data | | | |
|---|--|---------------------------------|-------------------------|
| article number | 13.32.05 | König ligature guide 5mm 19.5cm | |
| Unit | piece | Customs code (HS code) | 90189084 |
| GMDN | 45018 | GTIN | 04046826019686 |
| Risk category | 1r (EU MDR) | UMDNS-Code | 15-868 |
| Classification | 6 (EU MDR) | Base UDI-DI | 4046826;B;R001;03;ag;3J |
| Net weight (kg) | 0,036 | Operating instructions | GATDR00X |
| Latex free | yes | Preparation instructions | MPK03-1 |
|   | | | |
| Medicon media | 450.00.03 (general catalogue the art of surgery) 450.00.37 (catalogue ENT) | | |
| Complete picture | | Detail images | |
|    | | | |
| size | 5mm 19.5cm | Coating | |
| Storage conditions | | | |
| No special requirements | | | |
| Materials | | | |
| | Material number (DIN designation) | | |
| Main part Handle | 1.4305 (X8CrNiS18-9) 1.4301 (X5CrNi18-10) | | |