








## Product data sheet REF 13.32.07

11.12.2025

| General data  |   |                                 |                         |
|---|---|---------------------------------|-------------------------|
| article number  | 13.32.07  | König ligature guide 7mm 19.5cm |                         |
| Unit  | piece   | Customs code (HS code)          | 90189084                |
| GMDN  | 45018   | GTIN                            | 04046826019693          |
| Risk category   | 1r (EU MDR)   | UMDNS-Code                      | 15-868                  |
| Classification  | 6 (EU MDR)  | Base UDI-DI                     | 4046826;B;R001;03;ag;3J |
| Net weight (kg)   | 0,040   | Operating instructions          | GATDR00X                |
| Latex free  | yes   | Preparation instructions        | MPK03-1                 |
|     |   |                                 |                         |
| Medicon media   | 450.00.03 (general catalogue the art of surgery)<br>450.00.37 (catalogue ENT) |                                 |                         |
| Complete picture  |   | Detail images                   |                         |
|    |   |                                 |                         |
| size  | 7mm 19.5cm  | Coating                         |                         |
| Storage conditions  |   |                                 |                         |
| No special requirements   |   |                                 |                         |
| Materials   |   |                                 |                         |
|   | Material number (DIN designation)   |                                 |                         |
| Main part<br>Handle   | 1.4305 (X8CrNiS18-9)<br>1.4301 (X5CrNi18-10)                                  |                                 |                         |