



Product data sheet REF 14.44.67

25.08.2026

General data			
article number	14.44.67	Medicon HM wire twisting fcps 17cm	
Unit	piece	Customs code (HS code)	90189084
GMDN	32886	GTIN	04046826021023
Risk category	1r (EU MDR)	UMDNS-Code	15-221
Classification	6 (EU MDR)	Base UDI-DI	4046826;B;R002;02;ah;3V
Net weight (kg)	0,054	Operating instructions	GATDR00X
Latex free	yes	Preparation instructions	
 			
Medicon media	450.00.03 (general catalogue the art of surgery) 450.00.37 (catalogue ENT) 450.00.53 (catalogue cardio) 450.00.63 (catalogue Orthopaedics)		
Complete picture		Detail images	
			
size	17cm	Coating	Artificially created galvanic gold layer (max 0.04mm) on parts of the surface of instruments
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4021 (X20Cr13)		
Others	WSG_20-GO-55_CRTZ (RCoCr-C)		