



Product data sheet REF 15.04.35

03.06.2026

General data			
article number	15.04.35	Ti Bulldog clamp 3.5cm str V	
Unit	piece	Customs code (HS code)	90189084
GMDN	62470	GTIN	04046826022143
Risk category	1r (EU MDR)	UMDNS-Code	10-868
Classification	6 (EU MDR)	Base UDI-DI	4046826;B;R002;04;ah;4K
Net weight (kg)	0,006	Operating instructions	GATDR00X
Latex free	yes	Preparation instructions	MPK20-1



Medicon media	450.00.03 (general catalogue the art of surgery) 450.00.53 (catalogue cardio) 450.00.58 (catalogue neuro cranialis instr f cranial neurosurgery) 451.04.41 (brochure titanium instruments f Cardiac and Vascular Surgery)
---------------	---

Complete picture

Detail images



size	3.5cm	Coating	
------	-------	---------	--

Storage conditions

No special requirements

Materials

	Material number (DIN designation)
Main part	3.7035 (ASTM Grade 2)
Screw	1.4401 (X5CrNiMo17-12-2)
Others	1.4310 (X10CrNi18-8)