



Product data sheet REF 15.07.50

25.08.2026

General data			
article number	15.07.50	Ti Bulldog clamp 5.0cm strong cvd V	
Unit	piece	Customs code (HS code)	90189084
GMDN	62470	GTIN	04046826022358
Risk category	1r (EU MDR)	UMDNS-Code	10-868
Classification	6 (EU MDR)	Base UDI-DI	4046826;B;R002;04;ah;4K
Net weight (kg)	0,009	Operating instructions	GATDR00X
Latex free	yes	Preparation instructions	MPK20-1
			
Medicon media	450.00.03 (general catalogue the art of surgery) 450.00.53 (catalogue cardio) 450.00.58 (catalogue neuro cranialis instr f cranial neurosurgery) 451.04.41 (brochure titanium instruments f Cardiac and Vascular Surgery)		
Complete picture		Detail images	
			
size	5.0cm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	3.7035 (ASTM_Grade_2)		
Screw	1.4305 (X8CrNiS18-9)		
Others	1.4310 (X10CrNi18-8)		