



## Product data sheet REF 15.33.11

03.06.2026

| General data                                                                                                                                                                     |                                                  |                                        |                         |
|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------|----------------------------------------|-------------------------|
| article number                                                                                                                                                                   | 15.33.11                                         | Carrel (Reynolds) hemost fcps 11cm cvd |                         |
| Unit                                                                                                                                                                             | piece                                            | Customs code (HS code)                 | 90189084                |
| GMDN                                                                                                                                                                             | 62468                                            | GTIN                                   | 04046826022945          |
| Risk category                                                                                                                                                                    | 1r (EU MDR)                                      | UMDNS-Code                             | 10-865                  |
| Classification                                                                                                                                                                   | 6 (EU MDR)                                       | Base UDI-DI                            | 4046826;B;R002;04;ah;4K |
| Net weight (kg)                                                                                                                                                                  | 0,021                                            | Operating instructions                 | GATDR00X                |
| Latex free                                                                                                                                                                       | yes                                              | Preparation instructions               | MPK23-1                 |
| <div></div>    |                                                  |                                        |                         |
| Medicon media                                                                                                                                                                    | 450.00.03 (general catalogue the art of surgery) |                                        |                         |
| Complete picture                                                                                                                                                                 |                                                  | Detail images                          |                         |
| <div></div> |                                                  |                                        |                         |
| size                                                                                                                                                                             | 11cm                                             | Coating                                |                         |
| Storage conditions                                                                                                                                                               |                                                  |                                        |                         |
| No special requirements                                                                                                                                                          |                                                  |                                        |                         |
| Materials                                                                                                                                                                        |                                                  |                                        |                         |
|                                                                                                                                                                                  | Material number (DIN designation)                |                                        |                         |
| Main part                                                                                                                                                                        | 1.4021 (X20Cr13)                                 |                                        |                         |