



## Product data sheet REF 15.85.68

03.06.2026

| General data                                                                                                                                                                     |                                   |                                            |                         |
|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------|--------------------------------------------|-------------------------|
| article number                                                                                                                                                                   | 15.85.68                          | DeBaKey Bainbridge AT hemost fcps 18cm cvd |                         |
| Unit                                                                                                                                                                             | piece                             | Customs code (HS code)                     | 90189084                |
| GMDN                                                                                                                                                                             | 62468                             | GTIN                                       | 04046826023805          |
| Risk category                                                                                                                                                                    | 1r (EU MDR)                       | UMDNS-Code                                 | 10-865                  |
| Classification                                                                                                                                                                   | 6 (EU MDR)                        | Base UDI-DI                                | 4046826;B;R002;04;ah;4K |
| Net weight (kg)                                                                                                                                                                  | 0,035                             | Operating instructions                     | GATDR00X                |
| Latex free                                                                                                                                                                       | yes                               | Preparation instructions                   | MPK23-1                 |
| <div></div>    |                                   |                                            |                         |
| Medicon media                                                                                                                                                                    | 450.00.53 (catalogue cardio )     |                                            |                         |
| Complete picture                                                                                                                                                                 |                                   | Detail images                              |                         |
| <div></div> |                                   |                                            |                         |
| size                                                                                                                                                                             | 18cm                              | Coating                                    |                         |
| Storage conditions                                                                                                                                                               |                                   |                                            |                         |
| No special requirements                                                                                                                                                          |                                   |                                            |                         |
| Materials                                                                                                                                                                        |                                   |                                            |                         |
|                                                                                                                                                                                  | Material number (DIN designation) |                                            |                         |
| Main part                                                                                                                                                                        | 1.4021 (X20Cr13)                  |                                            |                         |