



## Product data sheet REF 16.15.92

11.12.2025

| General data                                                                                                                                                                  |                                                                                   |                                       |                         |
|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------|---------------------------------------|-------------------------|
| article number                                                                                                                                                                | 16.15.92                                                                          | Crafoord Sellors hemost fcps 24cm cvd |                         |
| Unit                                                                                                                                                                          | piece                                                                             | Customs code (HS code)                | 90189084                |
| GMDN                                                                                                                                                                          | 62468                                                                             | GTIN                                  | 04046826024475          |
| Risk category                                                                                                                                                                 | 1r (EU MDR)                                                                       | UMDNS-Code                            | 10-865                  |
| Classification                                                                                                                                                                | 6 (EU MDR)                                                                        | Base UDI-DI                           | 4046826;B;R002;04;ah;4K |
| Net weight (kg)                                                                                                                                                               | 0,082                                                                             | Operating instructions                | GATDR00X                |
| Latex free                                                                                                                                                                    | yes                                                                               | Preparation instructions              | MPK23-1                 |
| <div></div> |                                                                                   |                                       |                         |
| Medicon media                                                                                                                                                                 | 450.00.03 (general catalogue the art of surgery)<br>450.00.53 (catalogue cardio ) |                                       |                         |
| Complete picture                                                                                                                                                              |                                                                                   | Detail images                         |                         |
|          |                                                                                   |                                       |                         |
| size                                                                                                                                                                          | 24cm                                                                              | Coating                               |                         |
| Storage conditions                                                                                                                                                            |                                                                                   |                                       |                         |
| No special requirements                                                                                                                                                       |                                                                                   |                                       |                         |
| Materials                                                                                                                                                                     |                                                                                   |                                       |                         |
|                                                                                                                                                                               | Material number (DIN designation)                                                 |                                       |                         |
| Main part                                                                                                                                                                     | 1.4021 (X20Cr13)                                                                  |                                       |                         |