



## Product data sheet REF 16.20.18

25.08.2026

| General data                                                                                                                                                                  |                                                  |                                    |                         |
|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------|------------------------------------|-------------------------|
| article number                                                                                                                                                                | 16.20.18                                         | Pean hemost fcps short jaws 18.0cm |                         |
| Unit                                                                                                                                                                          | piece                                            | Customs code (HS code)             | 90189084                |
| GMDN                                                                                                                                                                          | 62468                                            | GTIN                               | 04046826024536          |
| Risk category                                                                                                                                                                 | 1r (EU MDR)                                      | UMDNS-Code                         | 10-865                  |
| Classification                                                                                                                                                                | 6 (EU MDR)                                       | Base UDI-DI                        | 4046826;B;R002;04;ah;4K |
| Net weight (kg)                                                                                                                                                               | 0,055                                            | Operating instructions             | GATDR00X                |
| Latex free                                                                                                                                                                    | yes                                              | Preparation instructions           | MPK23-1                 |
| <div></div> |                                                  |                                    |                         |
| Medicon media                                                                                                                                                                 | 450.00.03 (general catalogue the art of surgery) |                                    |                         |
| Complete picture                                                                                                                                                              |                                                  | Detail images                      |                         |
|          |                                                  |                                    |                         |
| size                                                                                                                                                                          | 18.0cm                                           | Coating                            |                         |
| Storage conditions                                                                                                                                                            |                                                  |                                    |                         |
| No special requirements                                                                                                                                                       |                                                  |                                    |                         |
| Materials                                                                                                                                                                     |                                                  |                                    |                         |
|                                                                                                                                                                               | Material number (DIN designation)                |                                    |                         |
| Main part<br>Screw                                                                                                                                                            | 1.4021 (X20Cr13)<br>1.4305 (X8CrNiS18-9)         |                                    |                         |