



## Product data sheet REF 16.33.23

25.08.2026

| General data    |             |                                      |                         |
|-----------------|-------------|--------------------------------------|-------------------------|
| article number  | 16.33.23    | Moynihan hemost fcps 23cm slight cvd |                         |
| Unit            | piece       | Customs code (HS code)               | 90189084                |
| GMDN            | 62468       | GTIN                                 | 04046826024796          |
| Risk category   | 1r (EU MDR) | UMDNS-Code                           | 10-865                  |
| Classification  | 6 (EU MDR)  | Base UDI-DI                          | 4046826;B;R002;04;ah;4K |
| Net weight (kg) | 0,065       | Operating instructions               | GATDR00X                |
| Latex free      | yes         | Preparation instructions             | MPK23-1                 |



Medicon media 450.00.03 (general catalogue the art of surgery)

## Complete picture

## Detail images



size 23cm Coating

## Storage conditions

No special requirements

## Materials

|           | Material number (DIN designation) |
|-----------|-----------------------------------|
| Main part | 1.4021 (X20Cr13)                  |