



## Product data sheet REF 16.35.43

25.08.2026

| General data                                                                                                                                                                                                                                                         |                                                  |                                    |                         |
|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------|------------------------------------|-------------------------|
| article number                                                                                                                                                                                                                                                       | 16.35.43                                         | Wertheim fcps hemost 1x2t 23cm cvd |                         |
| Unit                                                                                                                                                                                                                                                                 | piece                                            | Customs code (HS code)             | 90189084                |
| GMDN                                                                                                                                                                                                                                                                 | 62468                                            | GTIN                               | 04046826024901          |
| Risk category                                                                                                                                                                                                                                                        | 1r (EU MDR)                                      | UMDNS-Code                         | 10-865                  |
| Classification                                                                                                                                                                                                                                                       | 6 (EU MDR)                                       | Base UDI-DI                        | 4046826;B;R002;04;ah;4K |
| Net weight (kg)                                                                                                                                                                                                                                                      | 0,074                                            | Operating instructions             | GATDR00X                |
| Latex free                                                                                                                                                                                                                                                           | yes                                              | Preparation instructions           | MPK23-1                 |
| <div></div>                                                                                        |                                                  |                                    |                         |
| Medicon media                                                                                                                                                                                                                                                        | 450.00.03 (general catalogue the art of surgery) |                                    |                         |
| Complete picture                                                                                                                                                                                                                                                     |                                                  | Detail images                      |                         |
| <div></div> |                                                  |                                    |                         |
| size                                                                                                                                                                                                                                                                 | 23cm                                             | Coating                            |                         |
| Storage conditions                                                                                                                                                                                                                                                   |                                                  |                                    |                         |
| No special requirements                                                                                                                                                                                                                                              |                                                  |                                    |                         |
| Materials                                                                                                                                                                                                                                                            |                                                  |                                    |                         |
|                                                                                                                                                                                                                                                                      | Material number (DIN designation)                |                                    |                         |
| Main part                                                                                                                                                                                                                                                            | 1.4021 (X20Cr13)                                 |                                    |                         |
| Screw                                                                                                                                                                                                                                                                | 1.4305 (X8CrNiS18-9)                             |                                    |                         |