



## Product data sheet REF 16.53.20

11.12.2025

| General data                                                                                                                                                                                                                                                         |                                                                                   |                                 |                         |
|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------|---------------------------------|-------------------------|
| article number                                                                                                                                                                                                                                                       | 16.53.20                                                                          | Heiss hemost fcps 1x2t 21cm cvd |                         |
| Unit                                                                                                                                                                                                                                                                 | piece                                                                             | Customs code (HS code)          | 90189084                |
| GMDN                                                                                                                                                                                                                                                                 | 62468                                                                             | GTIN                            | 04046826025335          |
| Risk category                                                                                                                                                                                                                                                        | 1r (EU MDR)                                                                       | UMDNS-Code                      | 10-865                  |
| Classification                                                                                                                                                                                                                                                       | 6 (EU MDR)                                                                        | Base UDI-DI                     | 4046826;B;R002;04;ah;4K |
| Net weight (kg)                                                                                                                                                                                                                                                      | 0,047                                                                             | Operating instructions          | GATDR00X                |
| Latex free                                                                                                                                                                                                                                                           | yes                                                                               | Preparation instructions        | MPK23-1                 |
| <div></div>                                                                                        |                                                                                   |                                 |                         |
| Medicon media                                                                                                                                                                                                                                                        | 450.00.03 (general catalogue the art of surgery)<br>450.00.53 (catalogue cardio ) |                                 |                         |
| Complete picture                                                                                                                                                                                                                                                     |                                                                                   | Detail images                   |                         |
| <div></div> |                                                                                   |                                 |                         |
| size                                                                                                                                                                                                                                                                 | 21cm                                                                              | Coating                         |                         |
| Storage conditions                                                                                                                                                                                                                                                   |                                                                                   |                                 |                         |
| No special requirements                                                                                                                                                                                                                                              |                                                                                   |                                 |                         |
| Materials                                                                                                                                                                                                                                                            |                                                                                   |                                 |                         |
|                                                                                                                                                                                                                                                                      | Material number (DIN designation)                                                 |                                 |                         |
| Main part                                                                                                                                                                                                                                                            | 1.4021 (X20Cr13)                                                                  |                                 |                         |