



## Product data sheet REF 17.21.56

03.06.2026


| General data    |             |                                          |                         |
|-----------------|-------------|------------------------------------------|-------------------------|
| article number  | 17.21.56    | Schnidt (Sawtell) fcps 19cm strongly cvd |                         |
| Unit            | piece       | Customs code (HS code)                   | 90189084                |
| GMDN            | 62468       | GTIN                                     | 04046826026370          |
| Risk category   | 1r (EU MDR) | UMDNS-Code                               | 15-212                  |
| Classification  | 6 (EU MDR)  | Base UDI-DI                              | 4046826;B;R002;04;ah;4K |
| Net weight (kg) | 0,044       | Operating instructions                   | GATDR00X                |
| Latex free      | yes         | Preparation instructions                 | MPK23-1                 |



|               |                                                                                |
|---------------|--------------------------------------------------------------------------------|
| Medicon media | 450.00.03 (general catalogue the art of surgery)<br>450.00.37 (catalogue ENT ) |
|---------------|--------------------------------------------------------------------------------|

## Complete picture

## Detail images



|      |      |         |  |
|------|------|---------|--|
| size | 19cm | Coating |  |
|------|------|---------|--|

## Storage conditions

No special requirements

## Materials

|           | Material number (DIN designation) |
|-----------|-----------------------------------|
| Main part | 1.4021 (X20Cr13)                  |