



## Product data sheet REF 17.50.05

03.06.2026

| General data    |            |                          |                    |
|-----------------|------------|--------------------------|--------------------|
| article number  | 17.50.05   | towel clamp 5.5cm        |                    |
| Unit            | piece      | Customs code (HS code)   | 90189084           |
| GMDN            | 34961      | GTIN                     | 04046826026578     |
| Risk category   | 1 (EU MDR) | UMDNS-Code               | 10-902             |
| Classification  | 1 (EU MDR) | Base UDI-DI              | 4046826;N001;05;WQ |
| Net weight (kg) | 0,006      | Operating instructions   | GATDN00X           |
| Latex free      | yes        | Preparation instructions |                    |



|               |                                                                                |
|---------------|--------------------------------------------------------------------------------|
| Medicon media | 450.00.03 (general catalogue the art of surgery)<br>450.00.37 (catalogue ENT ) |
|---------------|--------------------------------------------------------------------------------|

## Complete picture

## Detail images



|      |       |         |  |
|------|-------|---------|--|
| size | 5.5cm | Coating |  |
|------|-------|---------|--|

## Storage conditions

No special requirements

## Materials

|           | Material number (DIN designation) |
|-----------|-----------------------------------|
| Main part | 1.4021 (X20Cr13)                  |