



Product data sheet REF 17.51.18

03.06.2026

General data			
article number	17.51.18	Doyen towel clamp 18cm	
Unit	piece	Customs code (HS code)	90189084
GMDN	34961	GTIN	04046826026615
Risk category	1r (EU MDR)	UMDNS-Code	10-902
Classification	6 (EU MDR)	Base UDI-DI	4046826;B;R002;04;ah;4K
Net weight (kg)	0,049	Operating instructions	GATDR00X
Latex free	yes	Preparation instructions	
			
Medicon media	450.00.03 (general catalogue the art of surgery)		
Complete picture		Detail images	
			
size	18cm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4021 (X20Cr13)		