



Product data sheet REF 17.54.10

03.06.2026

General data			
article number	17.54.10	Tohoku towel clamp 10cm	
Unit	piece	Customs code (HS code)	90189084
GMDN	34961	GTIN	04046826026622
Risk category	1 (EU MDR)	UMDNS-Code	10-902
Classification	1 (EU MDR)	Base UDI-DI	4046826;N001;05;WQ
Net weight (kg)	0,015	Operating instructions	GATDN00X
Latex free	yes	Preparation instructions	
			
Medicon media	450.00.03 (general catalogue the art of surgery)		
Complete picture		Detail images	
 			
size	10cm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4021 (X20Cr13)		