







Product data sheet REF 17.55.08

11.12.2025

General data			
article number	17.55.08	Backhaus towel clamp 8.5cm	
Unit	piece	Customs code (HS code)	90189084
GMDN	34961	GTIN	04046826026653
Risk category	1r (EU MDR)	UMDNS-Code	10-902
Classification	6 (EU MDR)	Base UDI-DI	4046826;B;R002;04;ah;4K
Net weight (kg)	0,017	Operating instructions	GATDR00X
Latex free	yes	Preparation instructions	MPK23-1
 			
Medicon media	450.00.03 (general catalogue the art of surgery) 450.00.37 (catalogue ENT) 450.00.40 (Katalog Nr 5 Dental) 450.00.62 (catalogue Medicon instruments Oral and Maxillofacial Surgery)		
Complete picture		Detail images	
			
size	8.5cm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4021 (X20Cr13)		