



## Product data sheet REF 17.72.20

03.06.2026

| General data    |            |                                         |                    |
|-----------------|------------|-----------------------------------------|--------------------|
| article number  | 17.72.20   | Presbyterian Hospital tubing clamp 20cm |                    |
| Unit            | piece      | Customs code (HS code)                  | 90189084           |
| GMDN            | 32146      | GTIN                                    | 04046826026851     |
| Risk category   | 1 (EU MDR) | UMDNS-Code                              | 10-875             |
| Classification  | 1 (EU MDR) | Base UDI-DI                             | 4046826;N001;05;WQ |
| Net weight (kg) | 0,079      | Operating instructions                  | GATDN00X           |
| Latex free      | yes        | Preparation instructions                |                    |



|               |                                                                                |
|---------------|--------------------------------------------------------------------------------|
| Medicon media | 450.00.03 (general catalogue the art of surgery)<br>450.00.37 (catalogue ENT ) |
|---------------|--------------------------------------------------------------------------------|

## Complete picture

## Detail images



|      |      |         |  |
|------|------|---------|--|
| size | 20cm | Coating |  |
|------|------|---------|--|

## Storage conditions

No special requirements

## Materials

|           | Material number (DIN designation) |
|-----------|-----------------------------------|
| Main part | 1.4021 (X20Cr13)                  |