



## Product data sheet REF 17.72.78

11.12.2025

| General data    |            |                                |                    |
|-----------------|------------|--------------------------------|--------------------|
| article number  | 17.72.78   | tubing clamp heavy patt 18.0cm |                    |
| Unit            | piece      | Customs code (HS code)         | 90189084           |
| GMDN            | 32146      | GTIN                           | 04046826026882     |
| Risk category   | 1 (EU MDR) | UMDNS-Code                     | 10-875             |
| Classification  | 1 (EU MDR) | Base UDI-DI                    | 4046826;N001;05;WQ |
| Net weight (kg) | 0,063      | Operating instructions         | GATDN00X           |
| Latex free      | yes        | Preparation instructions       |                    |



|               |                                                                                |
|---------------|--------------------------------------------------------------------------------|
| Medicon media | 450.00.03 (general catalogue the art of surgery)<br>450.00.37 (catalogue ENT ) |
|---------------|--------------------------------------------------------------------------------|

## Complete picture

## Detail images



|      |        |         |  |
|------|--------|---------|--|
| size | 18.0cm | Coating |  |
|------|--------|---------|--|

## Storage conditions

No special requirements

## Materials

|           | Material number (DIN designation) |
|-----------|-----------------------------------|
| Main part | 1.4021 (X20Cr13)                  |