




## Product data sheet REF 20.10.02

11.12.2025

| General data                                                                                                                                                                                                                                                |                                                                                                                          |                                |                         |
|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------|--------------------------------|-------------------------|
| article number                                                                                                                                                                                                                                              | 20.10.02                                                                                                                 | Meyering retr 18cm 5x16mm dent |                         |
| Unit                                                                                                                                                                                                                                                        | piece                                                                                                                    | Customs code (HS code)         | 90189084                |
| GMDN                                                                                                                                                                                                                                                        | 45918                                                                                                                    | GTIN                           | 04046826027735          |
| Risk category                                                                                                                                                                                                                                               | 1r (EU MDR)                                                                                                              | UMDNS-Code                     | 15-242                  |
| Classification                                                                                                                                                                                                                                              | 6 (EU MDR)                                                                                                               | Base UDI-DI                    | 4046826;B;R004;03;aj;5L |
| Net weight (kg)                                                                                                                                                                                                                                             | 0,037                                                                                                                    | Operating instructions         | GATDR00X                |
| Latex free                                                                                                                                                                                                                                                  | yes                                                                                                                      | Preparation instructions       | MPK03-1                 |
| <div></div>                                                                               |                                                                                                                          |                                |                         |
| Medicon media                                                                                                                                                                                                                                               | 450.00.03 (general catalogue the art of surgery)<br>450.00.37 (catalogue ENT )<br>450.00.60 (catalogue plastic surgery ) |                                |                         |
| Complete picture                                                                                                                                                                                                                                            |                                                                                                                          | Detail images                  |                         |
|    |                                                                                                                          |                                |                         |
| size                                                                                                                                                                                                                                                        | 18cm 5x16mm                                                                                                              | Coating                        |                         |
| Storage conditions                                                                                                                                                                                                                                          |                                                                                                                          |                                |                         |
| No special requirements                                                                                                                                                                                                                                     |                                                                                                                          |                                |                         |
| Materials                                                                                                                                                                                                                                                   |                                                                                                                          |                                |                         |
|                                                                                                                                                                                                                                                             | Material number (DIN designation)                                                                                        |                                |                         |
| Main part                                                                                                                                                                                                                                                   | 1.4021 (X20Cr13)                                                                                                         |                                |                         |