



## Product data sheet REF 20.12.03

11.12.2025

General data			
article number	20.12.03	Strandell Stille retr 19cm 6x18mm toothed insulated	
Unit	piece	Customs code (HS code)	90189084
GMDN	45918	GTIN	04046826539122
Risk category	1r (EU MDR)	UMDNS-Code	13-381
Classification	6 (EU MDR)	Base UDI-DI	4046826;B;R004;03;aj;5L
Net weight (kg)	0,000	Operating instructions	GATDR00X
Latex free	yes	Preparation instructions	MPK03-1
 			
Medicon media	453.03.01 (Blepharoplasty - Instruments for blepharoplasty Brochure)		
Complete picture		Detail images	
 			
size	19cm	Coating	powdered thermoplastic polyamide coating in various colors
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4021 (X20Cr13)		