



## Product data sheet REF 203.02.16

25.08.2026

| General data                                                                                                                                                                  |                                                                                                                                        |                               |                                                                                              |
|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------|-------------------------------|----------------------------------------------------------------------------------------------|
| article number                                                                                                                                                                | 203.02.16                                                                                                                              | SC surg scrs str 16.5cm bl/bl |                                                                                              |
| Unit                                                                                                                                                                          | piece                                                                                                                                  | Customs code (HS code)        | 90189084                                                                                     |
| GMDN                                                                                                                                                                          | 38727                                                                                                                                  | GTIN                          | 04046826171407                                                                               |
| Risk category                                                                                                                                                                 | 1r (EU MDR)                                                                                                                            | UMDNS-Code                    | 13-480                                                                                       |
| Classification                                                                                                                                                                | 6 (EU MDR)                                                                                                                             | Base UDI-DI                   | 4046826;B;R001;01;ag;2U                                                                      |
| Net weight (kg)                                                                                                                                                               | 0,063                                                                                                                                  | Operating instructions        | GATDR00X                                                                                     |
| Latex free                                                                                                                                                                    | yes                                                                                                                                    | Preparation instructions      |                                                                                              |
| <div></div> |                                                                                                                                        |                               |                                                                                              |
| Medicon media                                                                                                                                                                 | 450.00.58 (catalogue neuro cranialis instr f cranial neurosurgery)<br>450.00.59 (catalogue neuro spinalis instr f spinal neurosurgery) |                               |                                                                                              |
| Complete picture                                                                                                                                                              |                                                                                                                                        | Detail images                 |                                                                                              |
|                                                                                             |                                                                                                                                        |                               |                                                                                              |
| size                                                                                                                                                                          | 16.5cm                                                                                                                                 | Coating                       | Artificially created galvanic gold layer (max 0.04mm) on parts of the surface of instruments |
| Storage conditions                                                                                                                                                            |                                                                                                                                        |                               |                                                                                              |
| No special requirements                                                                                                                                                       |                                                                                                                                        |                               |                                                                                              |
| Materials                                                                                                                                                                     |                                                                                                                                        |                               |                                                                                              |
|                                                                                                                                                                               | Material number (DIN designation)                                                                                                      |                               |                                                                                              |
| Main part                                                                                                                                                                     | 1.4117 (X38CrMoV15)                                                                                                                    |                               |                                                                                              |
| Screw                                                                                                                                                                         | 1.4305 (X8CrNiS18-9)                                                                                                                   |                               |                                                                                              |