






Product data sheet REF 203.04.13

11.12.2025

| General data | | | |
|---|--|-------------------------------|--|
| article number | 203.04.13 | SC oper scrs str 13.0cm sh/sh | |
| Unit | piece | Customs code (HS code) | 90189084 |
| GMDN | 38727 | GTIN | 04046826171445 |
| Risk category | 1r (EU MDR) | UMDNS-Code | 13-480 |
| Classification | 6 (EU MDR) | Base UDI-DI | 4046826;B;R001;01;ag;2U |
| Net weight (kg) | 0,031 | Operating instructions | GATDR00X |
| Latex free | yes | Preparation instructions | MPK23-1 |
|   | | | |
| Medicon media | 450.00.03 (general catalogue the art of surgery) 450.00.37 (catalogue ENT) 451.03.23 (brochure medicut scrs) | | |
| Complete picture | | Detail images | |
|  | | | |
| size | 13.0cm | Coating | Artificially created galvanic gold layer (max 0.04mm) on parts of the surface of instruments |
| Storage conditions | | | |
| No special requirements | | | |
| Materials | | | |
| | Material number (DIN designation) | | |
| Main part Screw | 1.4117 (X38CrMoV15) 1.4305 (X8CrNiS18-9) | | |