



## Product data sheet REF 203.57.67

25.08.2026

| General data                                                                                                                                                                  |                                                                                                                                                                |                                   |                                                                                              |
|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------|----------------------------------------------------------------------------------------------|
| article number                                                                                                                                                                | 203.57.67                                                                                                                                                      | SC Plus Mayo Stille scrs cvd 17cm |                                                                                              |
| Unit                                                                                                                                                                          | piece                                                                                                                                                          | Customs code (HS code)            | 90189084                                                                                     |
| GMDN                                                                                                                                                                          | 38727                                                                                                                                                          | GTIN                              | 04046826172008                                                                               |
| Risk category                                                                                                                                                                 | 1r (EU MDR)                                                                                                                                                    | UMDNS-Code                        | 13-480                                                                                       |
| Classification                                                                                                                                                                | 6 (EU MDR)                                                                                                                                                     | Base UDI-DI                       | 4046826;B;R001;01;ag;2U                                                                      |
| Net weight (kg)                                                                                                                                                               | 0,059                                                                                                                                                          | Operating instructions            | GATDR00X                                                                                     |
| Latex free                                                                                                                                                                    | yes                                                                                                                                                            | Preparation instructions          | MPK23-1                                                                                      |
| <div></div> |                                                                                                                                                                |                                   |                                                                                              |
| Medicon media                                                                                                                                                                 | 450.00.03 (general catalogue the art of surgery)<br>450.00.37 (catalogue ENT )<br>450.00.60 (catalogue plastic surgery )<br>451.03.23 (brochure medicut scrs ) |                                   |                                                                                              |
| Complete picture                                                                                                                                                              |                                                                                                                                                                | Detail images                     |                                                                                              |
|                                                                                             |                                                                                                                                                                |                                   |                                                                                              |
| size                                                                                                                                                                          | 17cm                                                                                                                                                           | Coating                           | Artificially created galvanic gold layer (max 0.04mm) on parts of the surface of instruments |
| Storage conditions                                                                                                                                                            |                                                                                                                                                                |                                   |                                                                                              |
| No special requirements                                                                                                                                                       |                                                                                                                                                                |                                   |                                                                                              |
| Materials                                                                                                                                                                     |                                                                                                                                                                |                                   |                                                                                              |
|                                                                                                                                                                               | Material number (DIN designation)                                                                                                                              |                                   |                                                                                              |
| Main part                                                                                                                                                                     | 1.4117 (X38CrMoV15)                                                                                                                                            |                                   |                                                                                              |
| Screw                                                                                                                                                                         | 1.4305 (X8CrNiS18-9)                                                                                                                                           |                                   |                                                                                              |
| Others                                                                                                                                                                        | WSG_20-GO-55_CRTZ (RCoCr-C)                                                                                                                                    |                                   |                                                                                              |