



## Product data sheet REF 204.35.28

25.08.2026

| General data                                                                                                                                                                  |                                                                                                                                                                   |                                    |                                                                                              |
|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------------------------|----------------------------------------------------------------------------------------------|
| article number                                                                                                                                                                | 204.35.28                                                                                                                                                         | SC Metzenbaum slim scrs cvd 28.5cm |                                                                                              |
| Unit                                                                                                                                                                          | piece                                                                                                                                                             | Customs code (HS code)             | 90189084                                                                                     |
| GMDN                                                                                                                                                                          | 38727                                                                                                                                                             | GTIN                               | 04046826172671                                                                               |
| Risk category                                                                                                                                                                 | 1r (EU MDR)                                                                                                                                                       | UMDNS-Code                         | 13-492                                                                                       |
| Classification                                                                                                                                                                | 6 (EU MDR)                                                                                                                                                        | Base UDI-DI                        | 4046826;B;R001;01;ag;2U                                                                      |
| Net weight (kg)                                                                                                                                                               | 0,095                                                                                                                                                             | Operating instructions             | GATDR00X                                                                                     |
| Latex free                                                                                                                                                                    | yes                                                                                                                                                               | Preparation instructions           | MPK23-1                                                                                      |
| <div></div> |                                                                                                                                                                   |                                    |                                                                                              |
| Medicon media                                                                                                                                                                 | 450.00.03 (general catalogue the art of surgery)<br>450.00.53 (catalogue cardio )<br>450.00.60 (catalogue plastic surgery )<br>451.03.23 (brochure medicut scrs ) |                                    |                                                                                              |
| Complete picture                                                                                                                                                              |                                                                                                                                                                   | Detail images                      |                                                                                              |
|                                                                                             |                                                                                                                                                                   |                                    |                                                                                              |
| size                                                                                                                                                                          | 28.5cm                                                                                                                                                            | Coating                            | Artificially created galvanic gold layer (max 0.04mm) on parts of the surface of instruments |
| Storage conditions                                                                                                                                                            |                                                                                                                                                                   |                                    |                                                                                              |
| No special requirements                                                                                                                                                       |                                                                                                                                                                   |                                    |                                                                                              |
| Materials                                                                                                                                                                     |                                                                                                                                                                   |                                    |                                                                                              |
|                                                                                                                                                                               | Material number (DIN designation)                                                                                                                                 |                                    |                                                                                              |
| Main part                                                                                                                                                                     | 1.4117 (X38CrMoV15)                                                                                                                                               |                                    |                                                                                              |
| Screw                                                                                                                                                                         | 1.4305 (X8CrNiS18-9)                                                                                                                                              |                                    |                                                                                              |