



## Product data sheet REF 205.15.20

25.08.2026

| General data                                                                                                                                                                  |                                                                                                                                                                |                               |                                                                                              |
|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------|----------------------------------------------------------------------------------------------|
| article number                                                                                                                                                                | 205.15.20                                                                                                                                                      | SC micro 2000 scrs str 15.0cm |                                                                                              |
| Unit                                                                                                                                                                          | piece                                                                                                                                                          | Customs code (HS code)        | 90189084                                                                                     |
| GMDN                                                                                                                                                                          | 38727                                                                                                                                                          | GTIN                          | 04046826173258                                                                               |
| Risk category                                                                                                                                                                 | 1r (EU MDR)                                                                                                                                                    | UMDNS-Code                    | 15-621                                                                                       |
| Classification                                                                                                                                                                | 6 (EU MDR)                                                                                                                                                     | Base UDI-DI                   | 4046826;B;R001;01;ag;2U                                                                      |
| Net weight (kg)                                                                                                                                                               | 0,060                                                                                                                                                          | Operating instructions        | GATDR00X                                                                                     |
| Latex free                                                                                                                                                                    | yes                                                                                                                                                            | Preparation instructions      | MPK23-1                                                                                      |
| <div></div> |                                                                                                                                                                |                               |                                                                                              |
| Medicon media                                                                                                                                                                 | 450.00.03 (general catalogue the art of surgery)<br>450.00.37 (catalogue ENT )<br>450.00.60 (catalogue plastic surgery )<br>451.03.23 (brochure medicut scrs ) |                               |                                                                                              |
| Complete picture                                                                                                                                                              |                                                                                                                                                                | Detail images                 |                                                                                              |
|                                                                                             |                                                                                                                                                                |                               |                                                                                              |
| size                                                                                                                                                                          | 15.0cm                                                                                                                                                         | Coating                       | Artificially created galvanic gold layer (max 0.04mm) on parts of the surface of instruments |
| Storage conditions                                                                                                                                                            |                                                                                                                                                                |                               |                                                                                              |
| No special requirements                                                                                                                                                       |                                                                                                                                                                |                               |                                                                                              |
| Materials                                                                                                                                                                     |                                                                                                                                                                |                               |                                                                                              |
|                                                                                                                                                                               | Material number (DIN designation)                                                                                                                              |                               |                                                                                              |
| Main part                                                                                                                                                                     | 1.4021 (X20Cr13)                                                                                                                                               |                               |                                                                                              |
| Screw                                                                                                                                                                         | 1.4305 (X8CrNiS18-9)                                                                                                                                           |                               |                                                                                              |
| Others                                                                                                                                                                        | 1.4310 (X10CrNi18-8)                                                                                                                                           |                               |                                                                                              |