



## Product data sheet REF 214.36.11

11.12.2025

| General data                                                                                                                                                                  |                                                                                                                                                                |                                                     |                                                                                            |
|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------|--------------------------------------------------------------------------------------------|
| article number                                                                                                                                                                | 214.36.11                                                                                                                                                      | BI-Line SC Baby Metzenbaum scrs str 11.5cm Big Ring |                                                                                            |
| Unit                                                                                                                                                                          | piece                                                                                                                                                          | Customs code (HS code)                              | 90189084                                                                                   |
| GMDN                                                                                                                                                                          | 38727                                                                                                                                                          | GTIN                                                | 04046826174071                                                                             |
| Risk category                                                                                                                                                                 | 1r (EU MDR)                                                                                                                                                    | UMDNS-Code                                          | 13-492                                                                                     |
| Classification                                                                                                                                                                | 6 (EU MDR)                                                                                                                                                     | Base UDI-DI                                         | 4046826;B;R001;01;ag;2U                                                                    |
| Net weight (kg)                                                                                                                                                               | 0,025                                                                                                                                                          | Operating instructions                              | GATDR00X                                                                                   |
| Latex free                                                                                                                                                                    | yes                                                                                                                                                            | Preparation instructions                            | MPK23-1                                                                                    |
| <div></div> |                                                                                                                                                                |                                                     |                                                                                            |
| Medicon media                                                                                                                                                                 | 450.00.03 (general catalogue the art of surgery)<br>450.00.37 (catalogue ENT )<br>450.00.60 (catalogue plastic surgery )<br>451.03.23 (brochure medicut scrs ) |                                                     |                                                                                            |
| Complete picture                                                                                                                                                              |                                                                                                                                                                | Detail images                                       |                                                                                            |
|                                                                                            |                                                                                                                                                                |                                                     |                                                                                            |
| size                                                                                                                                                                          | 11.5cm                                                                                                                                                         | Coating                                             | Artificially created PVD layer of titanium aluminum carbon nitrite e.g. Magic Cut scissors |
|                                                                                                                                                                               |                                                                                                                                                                |                                                     |                                                                                            |
| Storage conditions                                                                                                                                                            |                                                                                                                                                                |                                                     |                                                                                            |
| No special requirements                                                                                                                                                       |                                                                                                                                                                |                                                     |                                                                                            |
| Materials                                                                                                                                                                     |                                                                                                                                                                |                                                     |                                                                                            |
|                                                                                                                                                                               | Material number (DIN designation)                                                                                                                              |                                                     |                                                                                            |
| Main part<br>Screw                                                                                                                                                            | 1.4117 (X38CrMoV15)                                                                                                                                            |                                                     |                                                                                            |
|                                                                                                                                                                               | 1.4305 (X8CrNiS18-9)                                                                                                                                           |                                                     |                                                                                            |