



## Product data sheet REF 214.75.23

03.06.2026

| General data                                                                                                                                                                  |                                                                                                                                                                                                 |                                     |                                                                                            |
|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------|--------------------------------------------------------------------------------------------|
| article number                                                                                                                                                                | 214.75.23                                                                                                                                                                                       | BI-Line SC Gorney Freeman scrs 23cm |                                                                                            |
| Unit                                                                                                                                                                          | piece                                                                                                                                                                                           | Customs code (HS code)              | 90189084                                                                                   |
| GMDN                                                                                                                                                                          | 38727                                                                                                                                                                                           | GTIN                                | 04046826174385                                                                             |
| Risk category                                                                                                                                                                 | 1r (EU MDR)                                                                                                                                                                                     | UMDNS-Code                          | 13-492                                                                                     |
| Classification                                                                                                                                                                | 6 (EU MDR)                                                                                                                                                                                      | Base UDI-DI                         | 4046826;B;R001;01;ag;2U                                                                    |
| Net weight (kg)                                                                                                                                                               | 0,070                                                                                                                                                                                           | Operating instructions              | GATDR00X                                                                                   |
| Latex free                                                                                                                                                                    | yes                                                                                                                                                                                             | Preparation instructions            | MPK23-1                                                                                    |
| <div></div> |                                                                                                                                                                                                 |                                     |                                                                                            |
| Medicon media                                                                                                                                                                 | 450.00.03 (general catalogue the art of surgery)<br>450.00.37 (catalogue ENT )<br>450.00.53 (catalogue cardio )<br>450.00.60 (catalogue plastic surgery )<br>451.03.23 (brochure medicut scrs ) |                                     |                                                                                            |
| Complete picture                                                                                                                                                              |                                                                                                                                                                                                 | Detail images                       |                                                                                            |
|                                                                                             |                                                                                                                                                                                                 |                                     |                                                                                            |
| size                                                                                                                                                                          | 23cm                                                                                                                                                                                            | Coating                             | Artificially created PVD layer of titanium aluminum carbon nitrite e.g. Magic Cut scissors |
| Storage conditions                                                                                                                                                            |                                                                                                                                                                                                 |                                     |                                                                                            |
| No special requirements                                                                                                                                                       |                                                                                                                                                                                                 |                                     |                                                                                            |
| Materials                                                                                                                                                                     |                                                                                                                                                                                                 |                                     |                                                                                            |
|                                                                                                                                                                               | Material number (DIN designation)                                                                                                                                                               |                                     |                                                                                            |
| Main part                                                                                                                                                                     | 1.4117 (X38CrMoV15)                                                                                                                                                                             |                                     |                                                                                            |
| Screw                                                                                                                                                                         | 1.4305 (X8CrNiS18-9)                                                                                                                                                                            |                                     |                                                                                            |