








Product data sheet REF 22.12.15

11.12.2025

| General data | | | |
|---|--|--------------------------|-------------------------|
| article number | 22.12.15 | Collin retr 16.5cm set/2 | |
| Unit | piece | Customs code (HS code) | 90189084 |
| GMDN | 45918 | GTIN | 04046826029265 |
| Risk category | 1r (EU MDR) | UMDNS-Code | 15-242 |
| Classification | 6 (EU MDR) | Base UDI-DI | 4046826;B;R004;03;aj;5L |
| Net weight (kg) | 0,092 | Operating instructions | GATDR00X |
| Latex free | yes | Preparation instructions | MPK12-1 |
|   | | | |
| Medicon media | 450.00.03 (general catalogue the art of surgery) 450.00.62 (catalogue Medicon instruments Oral and Maxillofacial Surgery) | | |
| Complete picture | | Detail images | |
|    | | | |
| size | 16.5cm | Coating | |
| Storage conditions | | | |
| No special requirements | | | |
| Materials | | | |
| | Material number (DIN designation) | | |
| Main part | 1.4301 (X5CrNi18-10) | | |