





Product data sheet REF 22.15.53

25.08.2026

General data			
article number	22.15.53	Sofield retr 20cm fig 3	
Unit	piece	Customs code (HS code)	90189084
GMDN	45918	GTIN	04046826029340
Risk category	1r (EU MDR)	UMDNS-Code	15-242
Classification	6 (EU MDR)	Base UDI-DI	4046826;B;R004;03;aj;5L
Net weight (kg)	0,071	Operating instructions	GATDR00X
Latex free	yes	Preparation instructions	MPK12-1
			
Medicon media	450.00.03 (general catalogue the art of surgery)		
Complete picture		Detail images	
  			
size	20cm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4021 (X20Cr13)		