



## Product data sheet REF 22.29.02

03.06.2026

| General data    |             |                             |                         |
|-----------------|-------------|-----------------------------|-------------------------|
| article number  | 22.29.02    | Volkman retr 2t bl 9mm 20cm |                         |
| Unit            | piece       | Customs code (HS code)      | 90189084                |
| GMDN            | 45918       | GTIN                        | 04046826029630          |
| Risk category   | 1r (EU MDR) | UMDNS-Code                  | 15-242                  |
| Classification  | 6 (EU MDR)  | Base UDI-DI                 | 4046826;B;R004;03;aj;5L |
| Net weight (kg) | 0,051       | Operating instructions      | GATDR00X                |
| Latex free      | yes         | Preparation instructions    | MPK12-1                 |



Medicon media 450.00.03 (general catalogue the art of surgery)

## Complete picture

## Detail images



size 9mm 20cm Coating

## Storage conditions

No special requirements

## Materials

|           | Material number (DIN designation) |
|-----------|-----------------------------------|
| Main part | 1.4021 (X20Cr13)                  |