



## Product data sheet REF 224.73.73

25.08.2026

| General data                                                                                                                                                                  |                                                                                                                                                                                                 |                          |                                                                                       |
|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------|---------------------------------------------------------------------------------------|
| article number                                                                                                                                                                | 224.73.73                                                                                                                                                                                       | MC Gorney scrs cvd 23cm  |                                                                                       |
| Unit                                                                                                                                                                          | piece                                                                                                                                                                                           | Customs code (HS code)   | 90189084                                                                              |
| GMDN                                                                                                                                                                          | 38727                                                                                                                                                                                           | GTIN                     | 04046826174972                                                                        |
| Risk category                                                                                                                                                                 | 1r (EU MDR)                                                                                                                                                                                     | UMDNS-Code               | 13-492                                                                                |
| Classification                                                                                                                                                                | 6 (EU MDR)                                                                                                                                                                                      | Base UDI-DI              | 4046826;B;R001;01;ag;2U                                                               |
| Net weight (kg)                                                                                                                                                               | 0,060                                                                                                                                                                                           | Operating instructions   | GATDR00X                                                                              |
| Latex free                                                                                                                                                                    | yes                                                                                                                                                                                             | Preparation instructions | MPK23-1                                                                               |
| <div></div> |                                                                                                                                                                                                 |                          |                                                                                       |
| Medicon media                                                                                                                                                                 | 450.00.03 (general catalogue the art of surgery)<br>450.00.37 (catalogue ENT )<br>450.00.53 (catalogue cardio )<br>450.00.60 (catalogue plastic surgery )<br>451.03.23 (brochure medicut scrs ) |                          |                                                                                       |
| Complete picture                                                                                                                                                              |                                                                                                                                                                                                 | Detail images            |                                                                                       |
|                                                                                             |                                                                                                                                                                                                 |                          |                                                                                       |
| size                                                                                                                                                                          | 23cm                                                                                                                                                                                            | Coating                  | Artificially created PVD layer of titanium aluminum nitrite and parts are gold plated |
| Storage conditions                                                                                                                                                            |                                                                                                                                                                                                 |                          |                                                                                       |
| No special requirements                                                                                                                                                       |                                                                                                                                                                                                 |                          |                                                                                       |
| Materials                                                                                                                                                                     |                                                                                                                                                                                                 |                          |                                                                                       |
|                                                                                                                                                                               | Material number (DIN designation)                                                                                                                                                               |                          |                                                                                       |
| Main part                                                                                                                                                                     | 1.4117 (X38CrMoV15)                                                                                                                                                                             |                          |                                                                                       |
| Screw                                                                                                                                                                         | 1.4305 (X8CrNiS18-9)                                                                                                                                                                            |                          |                                                                                       |
| Others                                                                                                                                                                        | WSG_20-GO-55_CRTZ (RCoCr-C)                                                                                                                                                                     |                          |                                                                                       |