



## Product data sheet REF 23.06.12

25.08.2026

| General data                                                                                                                                                                  |                                                                                                                                                           |                              |                         |
|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------|------------------------------|-------------------------|
| article number                                                                                                                                                                | 23.06.12                                                                                                                                                  | Langenbeck retr 14x33mm 21cm |                         |
| Unit                                                                                                                                                                          | piece                                                                                                                                                     | Customs code (HS code)       | 90189084                |
| GMDN                                                                                                                                                                          | 45918                                                                                                                                                     | GTIN                         | 04046826030377          |
| Risk category                                                                                                                                                                 | 1r (EU MDR)                                                                                                                                               | UMDNS-Code                   | 15-242                  |
| Classification                                                                                                                                                                | 6 (EU MDR)                                                                                                                                                | Base UDI-DI                  | 4046826;B;R004;03;aj;5L |
| Net weight (kg)                                                                                                                                                               | 0,080                                                                                                                                                     | Operating instructions       | GATDR00X                |
| Latex free                                                                                                                                                                    | yes                                                                                                                                                       | Preparation instructions     | MPK12-1                 |
| <div></div> |                                                                                                                                                           |                              |                         |
| Medicon media                                                                                                                                                                 | 450.00.03 (general catalogue the art of surgery)<br>450.00.37 (catalogue ENT )<br>450.00.53 (catalogue cardio )<br>450.00.60 (catalogue plastic surgery ) |                              |                         |
| Complete picture                                                                                                                                                              |                                                                                                                                                           | Detail images                |                         |
|                                                                                             |                                                                                                                                                           |                              |                         |
| size                                                                                                                                                                          | 14x33mm 21cm                                                                                                                                              | Coating                      |                         |
| Storage conditions                                                                                                                                                            |                                                                                                                                                           |                              |                         |
| No special requirements                                                                                                                                                       |                                                                                                                                                           |                              |                         |
| Materials                                                                                                                                                                     |                                                                                                                                                           |                              |                         |
|                                                                                                                                                                               | Material number (DIN designation)                                                                                                                         |                              |                         |
| Main part                                                                                                                                                                     | 1.4021 (X20Cr13)                                                                                                                                          |                              |                         |