



## Product data sheet REF 23.20.11

25.08.2026

| General data    |             |                          |                         |
|-----------------|-------------|--------------------------|-------------------------|
| article number  | 23.20.11    | Langenbeck retr 22x6mm   |                         |
| Unit            | piece       | Customs code (HS code)   | 90189084                |
| GMDN            | 45918       | GTIN                     | 04046826443139          |
| Risk category   | 1r (EU MDR) | UMDNS-Code               | 15-242                  |
| Classification  | 6 (EU MDR)  | Base UDI-DI              | 4046826;B;R004;03;aj;5L |
| Net weight (kg) | 0,094       | Operating instructions   | GATDR00X                |
| Latex free      | yes         | Preparation instructions | MPK12-1                 |



|               |                                                                          |
|---------------|--------------------------------------------------------------------------|
| Medicon media | 450.00.62 (catalogue Medicon instruments Oral and Maxillofacial Surgery) |
|---------------|--------------------------------------------------------------------------|

## Complete picture

## Detail images



|      |        |         |  |
|------|--------|---------|--|
| size | 22x6mm | Coating |  |
|------|--------|---------|--|

## Storage conditions

No special requirements

## Materials

|           |                                   |
|-----------|-----------------------------------|
|           | Material number (DIN designation) |
| Main part | 1.4021 (X20Cr13)                  |