



Product data sheet REF 23.25.09

25.08.2026

General data			
article number	23.25.09	Langenbeck Kocher retr 20x45mm 21cm	
Unit	piece	Customs code (HS code)	90189084
GMDN	45918	GTIN	04046826030698
Risk category	1r (EU MDR)	UMDNS-Code	15-242
Classification	6 (EU MDR)	Base UDI-DI	4046826;B;R004;03;aj;5L
Net weight (kg)	0,077	Operating instructions	GATDR00X
Latex free	yes	Preparation instructions	MPK12-1
			
Medicon media	450.00.03 (general catalogue the art of surgery) 450.00.37 (catalogue ENT) 450.00.46 (catalogue sinuscopy) 450.00.53 (catalogue cardio) 450.00.62 (catalogue Medicon instruments Oral and Maxillofacial Surgery)		
Complete picture		Detail images	
			
size	20x45mm 21cm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4021 (X20Cr13)		