



Product data sheet REF 23.37.20

25.08.2026

General data			
article number	23.37.20	Morris retr 50x55mm 22.5cm	
Unit	piece	Customs code (HS code)	90189084
GMDN	45918	GTIN	04046826031022
Risk category	1r (EU MDR)	UMDNS-Code	15-242
Classification	6 (EU MDR)	Base UDI-DI	4046826;B;R004;03;aj;5L
Net weight (kg)	0,147	Operating instructions	GATDR00X
Latex free	yes	Preparation instructions	MPK20-1



Medicon media 450.00.03 (general catalogue the art of surgery)

Complete picture

Detail images



size 50x55mm 22.5cm Coating

Storage conditions

No special requirements

Materials

	Material number (DIN designation)
Main part	1.4021 (X20Cr13)